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# MAINTENANCE AND SERVICE MANUAL

## Wheelchair Ramp – Residential Entrance

Modular System made of Steel / Aluminium

---

### 1. Purpose of this Manual

This manual describes the maintenance, inspection and care of modular wheelchair ramps for residential entrances made of steel and aluminium.

The purpose of this manual is to ensure the safety, functionality and long-term durability of the ramp throughout its entire service life.

This manual is intended for operators, property owners, property management companies, facility service providers and other responsible parties.

---

### 2. Product Overview

Modular wheelchair ramps for residential entrances consist of prefabricated modules that can be combined depending on the structural conditions on site. They can be permanently installed and are suitable for both indoor and outdoor use.

Typical components include:

- Ramp modules
- Platform(s)
- Guardrail and handrail systems, optional step modules
- Support and connection elements

The systems are available in steel as well as aluminium versions.

### 3. Intended Use

The ramp is intended exclusively for use by:

- Wheelchairs (manual and electric)
- Walking aids
- Pedestrians

Use with vehicles, machinery or for any other non-intended purpose is not permitted.

## 4. Operator Responsibilities

The operator is responsible for:

- Ensuring the safe condition of the ramp
- Performing regular visual and functional inspections
- Carrying out and documenting maintenance activities
- Immediate shutdown of the ramp in case of safety-relevant defects

## 5. Maintenance Intervals – Overview

Measure	Interval
General visual inspection	Monthly
Inspection of bolts & connections	Every 3–6 months
Cleaning of ramp surface	As required / seasonal
Winter operation (snow & ice)	Regularly
Inspection of handrails & guardrails	Monthly
Material condition / corrosion	At least annually

## **6. Regular Visual Inspection**

### **6.1 Ramp Surface**

- Free from dirt, leaves, ice or foreign objects
- Slip-resistant surface fully functional
- No sharp edges, cracks or deformation

### **6.2 Module Connections**

- All bolts and connection elements present
- No loose or damaged connections
- transitions between modules flush and stable

### **6.3 Support and Load-Bearing Elements**

- Stable, vertical positioning
- Secure contact with the ground
- No visible damage
- If step modules are installed, ensure secure fixing, stable positioning and intact tread surfaces.

## **7. Cleaning & Care (General)**

- Remove dirt and debris regularly
- Clean using water and mild detergent
- Do not use aggressive, acidic or chlorine-based cleaning agents

## **8. Material-Specific Instructions**

### **8.1 Wheelchair Ramp – Steel**

- Regularly inspect hot-dip galvanised surfaces
- Repair any damage to the galvanised coating professionally
- Do not use aggressive de-icing agents during winter

## 8.2 Wheelchair Ramp – Aluminium

- Keep aluminium surfaces clean to prevent abrasion
- Do not use abrasive or highly alkaline cleaning agents
- Regularly check all connection elements for tight fit

## 9. Handrails & Guardrails

- Secure fixing without movement
- No sharp edges or deformation
- Grip surfaces clean, hygienic and slip-resistant

**Damaged components must be repaired or replaced immediately.**

## 10. Safety Instructions

- The ramp must not be used if safety-relevant defects are present
- Repairs may only be carried out using suitable replacement parts
- Modifications or alterations may impair safety

## 11. Legal Notes

There is no legal requirement for mandatory TÜV inspection for wheelchair ramps of this type.

However, the operator is obliged to operate, inspect and maintain the ramp in accordance with applicable legal requirements and recognised rules of technology.

A voluntary inspection by an independent inspection body (e.g. TÜV) may be carried out, but is not mandatory.